


## Peeler - GUI

### Starting the program


- Open python console and import the peeler:  
`from Neuron_analysis_tool import peeler`
- This should open a matplotlib gui:



The left axis will show the voltage and the right axis will show the log voltage (for peeling).

### Loading a trace:

- Click on Load. and a file manager will open:



- Choose your trace. The trace need to be in the following format: 3/2 columns numpy readable, the first column is time, the second is current (you dont have to have this column), and the last one is


voltage. Like this:

```
peel.txt
~/PycharmProjects/Neuron_analysis_tool/Neuron_analysis_to...
Save

1 0.0000000000000000e+00 4.260253906250000000e+00 -7.275238037109375000e+01
2 2.0000000000000000e-02 4.150390625000000000e+00 -7.275327301025390625e+01
3 4.0000000000000000e-02 4.248046875000000000e+00 -7.275352478027343750e+01
4 5.9999999999999997e-02 4.199218750000000000e+00 -7.275377655029296875e+01
5 8.0000000000000000e-02 4.223632812500000000e+00 -7.275695800781250000e+01
6 1.0000000000000000e-01 4.211425781250000000e+00 -7.276039123535156250e+01
7 1.1999999999999999e-01 4.199218750000000000e+00 -7.276052093505859375e+01
8 1.4000000000000000e-01 4.174804687500000000e+00 -7.275848388671875000e+01
9 1.6000000000000000e-01 4.113769531250000000e+00 -7.275644683837890625e+01
10 1.7999999999999999e-01 4.138183593750000000e+00 -7.275835418701171875e+01
11 2.0000000000000000e-01 4.125976562500000000e+00 -7.275937652587890625e+01
12 2.2000000000000000e-01 4.248046875000000000e+00 -7.275873565673828125e+01
13 2.3999999999999999e-01 4.125976562500000000e+00 -7.275949859619140625e+01
14 2.6000000000000000e-01 4.174804687500000000e+00 -7.275785064697265625e+01
15 2.8000000000000000e-01 4.211425781250000000e+00 -7.275924682617187500e+01
16 2.9999999999999999e-01 4.064941406250000000e+00 -7.275886535644531250e+01
17 3.2000000000000000e-01 4.174804687500000000e+00 -7.275632476806640625e+01
18 3.4000000000000000e-01 4.260253906250000000e+00 -7.275924682617187500e+01
19 3.5999999999999998e-01 4.211425781250000000e+00 -7.275759124755859375e+01
20 3.8000000000000000e-01 4.260253906250000000e+00 -7.275453948974609375e+01
21 4.0000000000000000e-01 4.223632812500000000e+00 -7.275530242919921875e+01
22 4.1999999999999998e-01 4.223632812500000000e+00 -7.275314331054687500e+01
23 4.4000000000000000e-01 4.162597656250000000e+00 -7.275314331054687500e+01
24 4.6000000000000000e-01 4.211425781250000000e+00 -7.275225067138671875e+01
25 4.7999999999999998e-01 4.260253906250000000e+00 -7.275530242919921875e+01
26 5.0000000000000000e-01 4.138183593750000000e+00 -7.275365447998046875e+01
27 5.2000000000000000e-01 4.223632812500000000e+00 -7.275518035888671875e+01
28 5.4000000000000000e-01 4.199218750000000000e+00 -7.275428771972656250e+01
29 5.6000000000000000e-01 4.248046875000000000e+00 -7.275047302246093750e+01
30 5.7999999999999996e-01 4.138183593750000000e+00 -7.275377655029296875e+01
31 5.9999999999999997e-01 4.077148437500000000e+00 -7.275632476806640625e+01
32 6.1999999999999995e-01 4.211425781250000000e+00 -7.275556182861328125e+01
33 6.4000000000000000e-01 4.162597656250000000e+00 -7.274971008300781250e+01
34 6.6000000000000000e-01 4.211425781250000000e+00 -7.274971008300781250e+01
35 6.8000000000000000e-01 4.211425781250000000e+00 -7.275085449218750000e+01
36 7.0000000000000000e-01 4.150390625000000000e+00 -7.275339508056640625e+01
37 7.1999999999999997e-01 4.211425781250000000e+00 -7.274945831298828125e+01
38 7.3999999999999999e-01 4.223632812500000000e+00 -7.275009155273437500e+01
39 7.6000000000000000e-01 4.296875000000000000e+00 -7.275415802001953125e+01
40 7.8000000000000000e-01 4.162597656250000000e+00 -7.275327301025390625e+01
41 8.0000000000000000e-01 4.162597656250000000e+00 -7.275009155273437500e+01
42 8.2000000000000000e-01 4.199218750000000000e+00 -7.275327301025390625e+01
43 8.3999999999999996e-01 4.174804687500000000e+00 -7.275848388671875000e+01
44 8.5999999999999998e-01 4.284667968750000000e+00 -7.275848388671875000e+01
45 8.8000000000000000e-01 4.187011718750000000e+00 -7.275657653808593750e+01
46 9.0000000000000000e-01 4.162597656250000000e+00 -7.276039123535156250e+01
47 9.2000000000000000e-01 4.260253906250000000e+00 -7.276026153564453125e+01
48 9.4000000000000000e-01 4.125976562500000000e+00 -7.275644683837890625e+01
49 9.5999999999999996e-01 4.150390625000000000e+00 -7.275746917724609375e+01
50 9.7999999999999998e-01 4.211425781250000000e+00 -7.275606536865234375e+01
51 1.0000000000000000e+00 4.235839843750000000e+00 -7.276026153564453125e+01
52 1.0200000000000000e+00 4.174804687500000000e+00 -7.275988006591796875e+01
53 1.0400000000000000e+00 4.260253906250000000e+00 -7.275733947753906250e+01

Plain Text Tab Width: 8 Ln 1, Col 76 INS
```

- Click Open and the trace will load into view:



### Pre peeling setting

- Before peeling exponents, you can use the following options:
  - Cut - cut your point of interest from the trace (After choosing cut you need to choose 2 points on the trace for start and end).
  - Flip - flip the y axis (peeling is only acceptable for positive values so if you used negative current injection you need to Flip so that the voltage response will be positive).
  - Down/Up - shift all the trace up or down by a constant factor (we set this to 0.05 mV) but you can change it in peeler.py script SHIFT\_FACTOR.


**Notice that shifting the voltage can affect your results (see).**

- Return - go back to the previous step.


Example:

For the top screenshot:


Flip:



Cut:




Down:



### Peeling

- When you finished presenting the trace parameter you can star peeling.
- Choose the peel option and choose a linear part in the log-scale axis (the right plot).


- After choosing the 2 points you will get the following plot:



The red dots show your selected points and the line shows the line fit in log scale. The extracted time constant ( $\tau$ ) will appear at the plot title.

- To confirm click OK, to change selection click Return.
- After clicking OK 3 things will happen:
  - a. The peeled voltage will be reduced from the original voltage on both the voltage and log voltage axis.
  - b. The pre peeling setting will disappear.
  - c. The time constant you peeled so far will be listed on the voltage plot.






- You can continue to peel the nex exponent by clicking peel again.

### Saving the result

- You can click on Save that will open the file manager. Choose the location where you want to save the file and click save.

- The plot of the fitting:



Will be saved to the location you chose, and also a pickle object that you can load via Open saved button.

- The plots to the left show the cutted voltage part (solid black) and the peeling result in dashed red (how much of the signal is explained).
- The plot to the right shows the log voltage in (solid black) and the peeling lines in colors.

For any question please send an email to [yonil.leibner@mail.huji.ac.il](mailto:yonil.leibner@mail.huji.ac.il)